


ALOVA® 100



Braden scale score



KEY FEATURES

- ▶ Single-piece mattress made up of:
 - 6 cm molded visco elastic foam top layer: density $\geq 75 \text{ Kg/m}^3$
 - 8 cm High Resilient (HR) foam base: density $\geq 34 \text{ Kg/m}^3$
- ▶ Patient weight: 30 to 150 Kg
- ▶ Dimensions: 195 x 97 x 14 cm (76.7 x 38.1 x 5.5 inches)
- ▶ Mattress weight: 15 Kg
- ▶ Cover available: Promust CIC

KEY BENEFITS

- ▶ Visco elastic foam provides excellent pressure re-distribution properties: homogeneous immersion of the body into the foam, thus increasing the bearing surface and reducing pressure applied in areas at risk of pressure ulcers.
- ▶ 100 cm width for larger sleeping space and comfort.
- ▶ Visco elastic foam provides patient comfort.

INDICATIONS

HELP IN PREVENTING PRESSURE ULCERS FOR PATIENT AT MEDIUM TO HIGH RISK OF PRESSURE ULCER DEVELOPMENT, BEDRIDDEN MORE THAN 15 HOURS A DAY.

HELP IN TREATING EXISTING PRESSURE ULCERS FROM STADE 1 TO 2 COMBINED WITH TECHNICAL AID DEVICES.

PATIENT WEIGHT: 30 TO 150 KG/ 4.7 TO 23.6 ST

- ▶ Single-piece mattress made up of visco elastic foam top layer on High Resilient (HR) polyurethane foam
- ▶ Removable and washable integral cover

DESCRIPTION	ITEM REFERENCES ACCORDING TO PACKAGING		FIRE STANDARDS
	CARDBOX	PALLET BOX	
ALOVA® 100 Promust CIC cover	VMA100-CIC	VMA100-CIC-BOX	EN 597-1 and 2 GPEM D1-90 classe C

3 years warranty

Cleaning instructions: Promust CIC cover



FP000064/01.2018

Read carefully the product instructions before use

Photos not contractually binding



WINCARE Group - 200 rue Charles Tellier - Actiparc de Grézan - 30034 Nîmes Cedex 1 France
Ph : +33 (0)4 66 02 15 15 - Fax : +33 (0)4 66 02 15 00 - Email : contact.askle@winnicare.fr - www.winnicare.fr
In the interest of constantly improving its products, the WINCARE Group reserves the rights to modify the technical specifications of any product without prior notice
Class I medical device in accordance with European Directive 93/42/EEC

