

Luna hand control and charger

IPx7 screw lock connector:
A secure and reliable mechanism to hook up the hand control to the Luna lift.

- Luna hand control has ergonomic design and exceptional functionality.
- Luna charger can be mounted on every wall. This ensures a permanent and stable charging station.
- Luna charger can also be mounted on other Ergolet products, such as Luna Trolley or Mobile Flex 2- and 4-gantry systems.

Permanent and stable charging

Scan QR and see Luna video, product sheet, quickguide or user manual



Hand control models

- The hand control is now available in IPx4 and IPx7 protection class. IPx4 tolerates exposure to showering. IPx7 tolerates full immersion in water.
- The IPx4 hand control has optional top or bottom cable exit. The IPx7 hand control has bottom cable exit.



Model (class): IPx4
Cable exit: BOTTOM EXIT

Cable length: 1320mm (item no. 015-03101)
Cable length: 1620mm (item no. 015-03102)
Cable length: 1900mm (item no. 015-03103)

Tolerates: Exposure to showering

Model (class): IPx7
Cable exit: BOTTOM EXIT

Cable length: 1320mm (015-03109)

Tolerates: Full immersion in water

Model (class): IPx4, 2 channels
Cable exit: BOTTOM EXIT

Cable length: 1320mm (item no. 015-03106)
Cable length: 1620mm (item no. 015-03107)
Cable length: 1900mm (item no. 015-03108)

Tolerates: Exposure to showering

Model (class): IPx4
Cable exit: TOP EXIT

Cable length: 1320mm (item no. 015-03105)

Tolerate: Exposure to showering

Model (class): IPx4, 2 channels
Cable exit: TOP EXIT

Cable length: 1320mm (item no. 015-03110)
Cable length: 1620mm (item no. 015-03111)
Cable length: 1900mm (item no. 015-03112)

Tolerates: Exposure to showering



2 channels:
 additional buttons
 for additional
 equipment, such
 as Ergo Trainer or
 electric spreader bar

All protection class tests are carried out by TÜV.



MDD93/42/EEC
 Class I

