# Flex-Lift<sup>™</sup>/130

A mobile floor lift you can carry



### Primary use:

- Lift users, who travel a lot, either job-related or on vacation
- As a second lift for disabled children, who already use a ceiling lift
- For community nurses on house visits

With its excellent lifting range the Flex-Lift allows patients to be lifted from the floor.

This ergonomic mobile floor lift can be folded in seconds for easy and discreet storage.

It fits easily into the trunk of a standard car and can be carried as two separate bags weighing only 16 and 15 kg.

## Product Highlights:

- Max. user weight 130 kg. (286 lbs / 20,66 st.)
- Electrical lifting.
- Lifts from the floor. Lifting range of 116 cm (46")
- The chassis folds flat for convenient storing
- Can be assembled and disassembled in minutes, without the use of tools
- Can be carried effortlessly by one person
- Fits into a standard size trunk making for versatile use in many situations









Good lifting height.



Foldable and easy to store. Takes up little space.



Carrying bags with wheels Transport Bags: protects the lift and makes it easier to carry.



Fits into the trunk of a standard car.

| Lift                      | Order No.           |
|---------------------------|---------------------|
| Flex-Lift                 | 010-02001UK         |
| Accessories               |                     |
| Carrying bags with wheels | 017-02001           |
| Slings                    | Consult your dealer |
| Extra battery pack        | 018-02300           |

## **TECHNICAL SPECIFICATION**

### Flex-Lift, assembled:

| Height:                         | 117 - 174 cm | (46-69")              |
|---------------------------------|--------------|-----------------------|
| Hoisting range:                 | 116 cm       | (46")                 |
| Spreader bar, lowest position:  | 46 cm        | (18")                 |
| Spreader bar, highest position: | 162 cm       | (64")                 |
| Width, folded:                  | 28 cm        | (11")                 |
| Width, outside:                 | 77,5 cm      | (26")                 |
| Width, inside:                  | 66 cm        | (26")                 |
| Height of undercarriage:        | 11,5 cm      | (4,5")                |
| Clearance, undercarriage:       | 5,5 cm       | (2,2")                |
| Total weight:                   | 23,9 kg      | (52,6 lbs.)           |
| Max user weight:                | 130 kg       | (286 lbs. / 20,66 st) |

Battery:

Battery 24 volt

Battery capacity: 25 – 30 lifting cycles

Charging: Through hand control at external charging station

Charging time: 4-5 hours

Undercarriage:

 Length:
 113 cm
 (44,5")

 Width, folded:
 22,5 cm
 (9")

 Weight:
 12,6 kg
 (27,7 lbs.)

 Measurements in bag:
 128x32x34 cm
 46"x12"x13"

 Weight in bag:
 16,8 Kg
 (37,0 lbs.)

Mast:

 Length:
 110 cm
 (43")

 Weight:
 11,3 kg
 (24,9 lbs.)

 Measurements in bag:
 128x32x34 cm
 46"x12"x13"

 Weight in bag:
 15,5 Kg
 (34,1 lbs.)

Flex-Lift has been tested and approved according to EN 10535, EN 60601-1, EN 60601-1-2 Flex-Lift carries the CE-mark

Flex-Lift is patented - EP 1383455

Consult the user manual for further specifications. All figures and measurements are for guidance only. Winncare Nordic ApS reserves the right to changes without notice.

Product photos may include accessories. Consult your dealer.

Ergolet The Courtyard Wellington Road North, Stockport Cheshire SK4 1HT Tel: 0161 477 7900 Fax: 0161 480 7733 e-mail: info@ergolet.co.uk

